

To analyze the tubulin dynamics during spindle formation, *Drosophila* S2 cells were transfected with mEos2-labeled tubulin. By local irradiation with UV laser light (405nm) the green fluorescent mEos2 photoconverts to a red fluorescent dye. In contrast to photobleaching, where the fluorescent signal is locally depleted, the photoswitching allowed the independent observation of two different tubulin fractions in space and time.

[https://rapp-opto.com/wp-content/uploads/2018/10/thorn-ucsf-spindle-elongation\\_x264.mp4](https://rapp-opto.com/wp-content/uploads/2018/10/thorn-ucsf-spindle-elongation_x264.mp4)

Movie1: Photoconversion of mEos2-labeled tubulin in the spindle of a *Drosophila* S2 cell. The video is sped up 20-fold from real time.

**Setup:**

- Microscope: Standard widefield microscope
- Objective: 100x NA 1.4
- 405 & 473 nm diode laser

**Rapp OptoElectronic Components:**

- UGA-40 – point scanning device (integrated in  $\mu$ -manager)

**Data taken from:**

Kurt's Microscopy Blog

<http://nic.ucsf.edu/blog/2014/04/photobleaching-and-photoactivation/>

Kurt Thorn (1) & Nico Stuurman (2)

(1) Nikon Imaging Center (NIC) at University of California - San Francisco (UCSF)

(2) Vale Lab at University of California - San Francisco (UCSF)